Print this Page for Your Records

Close Window

** SESSION **

Neural Control and Balance

Presentation Number: 2439

Title:

The Effects of 4 Different Acute Whole Body Vibration Exposures upon Indicies

of Counter Movement Vertical Jump Performance

Presentation Start:

6/2/2006 8:30:00 AM

Presentation End:

6/2/2006 9:30:00 AM

Topical Category:

405. neural control of movement and balance

Authors:

Hugh S. Lamont, Mike G. Bemben, FACSM, Joel Cramer, Adrien Gayaud, Luke S.

Acree. University of Oklahoma, Norman, OK.

Email: hslamont@ou.edu

The use of whole body low frequency vibration to acutely excite the neuromuscular system is gaining popularity. PURPOSE: To investigate the affects of four different whole body vibration applications upon indices of vertical jump performance. Methods: Twenty-one recreationally trained males (179.1 ± 5.6 cm in height, 84.7 ± 11.4 kg in weight and 26 ± 12 vrs of age) participated in two, non-consecutive visits. Subjects performed 12 counter movement vertical jumps (CMVJ) upon a Just Jump System switch mat in a standardized position. The CMVJ'S were performed in four series [pre-vibration (T1), two minute post-vibration (T2), five-minutes and 30 seconds post-vibration (T3), 11 minutes postvibration (T4)] of three trials. Vibration treatments were given at a frequency of either 30 Hz or 50 Hz (continuous or intermittent) for a total of 30 seconds, for a total of four conditions. On each visit, subjects received two randomly assigned vibration treatments. Results: Jump height after intermittent vibration of 50 Hz was greater than after continuous vibration of 30 Hz (p = 0.039). Overall, differences were found among the four conditions and the four series for peak power (p = 0.022 and p = 0.035, respectively), relative peak power (p = 0.028 and p = 0.036, respectively), FITRO® dyne mean power (p = 0.012 and p = 0.003, respectively), and Lewis nomogram mean power (p = 0.014 and p = 0.018, respectively). The two-minute post-vibration trial had significantly higher mean power than all other trials (p < 0.05). Additionally, the pre-vibration trial yielded lower FITRO® dyne mean power values than the trials post-vibration (p = 0.035). Lastly, the intermittent vibration at 50 Hz yielded greater mean power than the continuous bouts of vibration (p < 0.05). Conclusion: The use of an intermittent vibration protocol at a frequency of 50 Hz appears to be more effective than continuous vibration at a lower frequency of 30 Hz at enhancing acute indices of vertical jump performance.

OASIS - Online Abstract Submission and Invitation System™ ©1996-2006, Coe-Truman Technologies, Inc.